

28 Liter Stainless Steel Open Tank



T28S

Key Specifications

Description	Stainless Steel Open Tank, 28L
Working Temperature Range °C	Up to 150°
Working Access (L x W x D) (cm)	44.5 x 24.1 x 20.3

Key Features

IMMERSION CIRCULATOR & OPEN TANKS

Key Features

- Large, universal icon and English display
- On-screen prompts
- Single-point calibration capability

Common Specifications

Working Temperature Range °F	Up to 302°
Working Temperature Range °C	Up to 150°
Reservoir Capacity (gallons)	7.4

Reservoir Capacity (liters)	28
Reservoir/Tank Material	Stainless Steel
Working Access (L x W x D) (inches)	17.5 x 9.5 x 8
Working Access (L x W x D) (cm)	44.5 x 24.1 x 20.3
Overall Dimensions (L x W x H) (inches)	20.87 x 12.87 x 9.62 in
Overall Dimensions (L x W x H) (cm)	53 x 32.7 x 24.4 cm
Catalog Page Number	67

60 Hz Only

Part Number T28SA1

50 Hz Only



© 2015-2020 **PolyScience** 6600 W. Touhy Avenue, Niles, Illinois 60714-4516 USA

Phone: +1(847) 647-0611 or +1(800) 229-7569 **Fax:** +1(847) 647-1155

*/

Source URL (modified on 12/27/2019 - 12:14): <https://polyscience.com/circulating-baths/28-liter-stainless-steel-open-tank>