

8 Liter Polycarbonate Open Tank



T08P

Key Specifications

Description	Polycarbonate Open Tank, 8L
Working Temperature Range °C	Up to 85°
Working Access (L x W x D) (cm)	28.4 x 15.6 x 20.3

Key Features

IMMERSION CIRCULATOR & OPEN TANKS

Key Features

- Large, universal icon and English display
- On-screen prompts
- Single-point calibration capability

Common Specifications

Working Temperature Range °F	Up to 185°
Working Temperature Range °C	Up to 85°
Reservoir Capacity (gallons)	2.1

Reservoir Capacity (liters)	8
Reservoir/Tank Material	Polycarbonate
Working Access (L x W x D) (inches)	11.2 x 6.13 x 8
Working Access (L x W x D) (cm)	28.4 x 15.6 x 20.3
Overall Dimensions (L x W x H) (inches)	11.9 x 6.87 x 9.13 in
Overall Dimensions (L x W x H) (cm)	30.2 x 17.4 x 23.2 cm
Catalog Page Number	67

60 Hz Only

Part Number T08PA1

50 Hz Only



© 2015-2020 **PolyScience** 6600 W. Touhy Avenue, Niles, Illinois 60714-4516 USA

Phone: +1(847) 647-0611 or +1(800) 229-7569 **Fax:** +1(847) 647-1155

*/

Source URL (modified on 09/23/2020 - 16:12): <https://polyscience.com/circulating-baths/t08p-8-liter-polycarbonate-open-tank>