

11 Liter Polycarbonate Open Tank



T11P

Key Specifications

| | |
|--|------------------------------|
| Description | Polycarbonate Open Tank, 11L |
| Working Temperature Range °C | Up to 85° |
| Working Access (L x W x D) (cm) | 38.9 x 15.6 x 20.3 |

Key Features

IMMERSION CIRCULATOR & OPEN TANKS

Key Features

- Large, universal icon and English display
- On-screen prompts
- Single-point calibration capability

Common Specifications

| | |
|-------------------------------------|------------|
| Working Temperature Range °F | Up to 185° |
| Working Temperature Range °C | Up to 85° |
| Reservoir Capacity (gallons) | 2.9 |

| | |
|--|-----------------------|
| Reservoir Capacity (liters) | 11 |
| Reservoir/Tank Material | Polycarbonate |
| Working Access (L x W x D) (inches) | 15.3 x 6.13 x 8 |
| Working Access (L x W x D) (cm) | 38.9 x 15.6 x 20.3 |
| Overall Dimensions (L x W x H) (inches) | 16 x 6.87 x 9.13 in |
| Overall Dimensions (L x W x H) (cm) | 40.6 x 17.4 x 23.2 cm |
| Catalog Page Number | 67 |

60 Hz Only

Part Number T11PA1

50 Hz Only



© 2015-2020 **PolyScience** 6600 W. Touhy Avenue, Niles, Illinois 60714-4516 USA

Phone: +1(847) 647-0611 or +1(800) 229-7569 **Fax:** +1(847) 647-1155

*/

Source URL (modified on 08/05/2020 - 13:00): <https://polyscience.com/circulating-baths/t11p-11-liter-polycarbonate-open-tank>