

# 14 Liter Polycarbonate Open Tank



T14P

## Key Specifications

<b>Description</b>	Polycarbonate Open Tank, 14L
<b>Working Temperature Range °C</b>	Up to 85°
<b>Working Access (L x W x D) (cm)</b>	49.3 x 15.6 x 20.3

## Key Features

### IMMERSION CIRCULATOR & OPEN TANKS

#### Key Features

- Large, universal icon and English display
- On-screen prompts
- Single-point calibration capability

## Common Specifications

<b>Working Temperature Range °F</b>	Up to 185°
<b>Working Temperature Range °C</b>	Up to 85°
<b>Reservoir Capacity (gallons)</b>	3.7

<b>Reservoir Capacity (liters)</b>	14
<b>Reservoir/Tank Material</b>	Polycarbonate
<b>Working Access (L x W x D) (inches)</b>	19.4 x 6.13 x 8
<b>Working Access (L x W x D) (cm)</b>	49.3 x 15.6 x 20.3
<b>Overall Dimensions (L x W x H) (inches)</b>	20.2 x 6.87 x 9.13 in
<b>Overall Dimensions (L x W x H) (cm)</b>	51.3 x 17.4 x 23.2 cm
<b>Catalog Page Number</b>	67

## 60 Hz Only

**Part Number** T14PA1

## 50 Hz Only



© 2015-2020 **PolyScience** 6600 W. Touhy Avenue, Niles, Illinois 60714-4516 USA

**Phone:** +1(847) 647-0611 or +1(800) 229-7569 **Fax:** +1(847) 647-1155

\*/

**Source URL (modified on 08/05/2020 - 16:10):** <https://polyscience.com/circulating-baths/t14p-14-liter-polycarbonate-open-tank>