

## 23 Liter Polycarbonate Open Tank



T23P

### Key Specifications

Description	Polycarbonate Open Tank, 23L
Working Temperature Range °C	Up to 85°
Working Access (L x W x D) (cm)	41.7 x 30.5 x 20.3

### Key Features

#### IMMERSION CIRCULATOR & OPEN TANKS

##### Key Features

- Large, universal icon and English display
- On-screen prompts
- Single-point calibration capability

### Common Specifications

Working Temperature Range °F	Up to 185°
Working Temperature Range °C	Up to 85°
Reservoir Capacity (gallons)	6.1

<b>Reservoir Capacity (liters)</b>	23
<b>Reservoir/Tank Material</b>	Polycarbonate
<b>Working Access (L x W x D) (inches)</b>	16.4 x 12 x 8
<b>Working Access (L x W x D) (cm)</b>	41.7 x 30.5 x 20.3
<b>Overall Dimensions (L x W x H) (inches)</b>	17.9 x 12.75 x 9.13 in
<b>Overall Dimensions (L x W x H) (cm)</b>	45.5 x 32.4 x 23.2 cm
<b>Catalog Page Number</b>	67

## 60 Hz Only

Part Number T23PA1

## 50 Hz Only



© 2015-2020 **PolyScience** 6600 W. Touhy Avenue, Niles, Illinois 60714-4516 USA

**Phone:** +1(847) 647-0611 or +1(800) 229-7569 **Fax:** +1(847) 647-1155

\*/

**Source URL (modified on 08/05/2020 - 16:39):** <https://polyscience.com/circulating-baths/t23p-23-liter-polycarbonate-open-tank>